

Authors:

Alberto Esquenazi, MD Mukul Talaty, PhD Andrew Packel, PT, NCS Michael Saulino, MD, PhD

Affiliations:

From MossRehab, Elkins Park, Pennsylvania.

Correspondence:

All correspondence and requests for reprints should be addressed to: Alberto Esquenazi, MD, 60 Township Line Road, Gait & Motion Analysis Laboratory, Elkins Park, PA 19027.

Disclosures:

Supported, in part, by Argo Medica and MossRehab. The ReWalk devices used were provided by the study sponsor. The device manufacturer did not have any role in the subject selection, access to the data and was not involved in the data analysis or publication.

Presented at the 2010 Annual Assembly of the American Academy of Physical Medicine and Rehabilitation in Seattle, WA (preliminary version).

Financial disclosure statements have been obtained, and no conflicts of interest have been reported by the authors or by any individuals in control of the content of this article.

0894-9115/12/0000-0000/0

American Journal of Physical Medicine & Rehabilitation Copyright * 2012 by Lippincott Williams & Wilkins

DOI: 10.1097/PHM.0b013e318269d9a3

SCI

The ReWalk Powered Exoskeleton to Restore Ambulatory Function to Individuals with Thoracic-Level Motor-Complete Spinal Cord Injury

ABSTRACT

Esquenazi A, Talaty M, Packel A, Saulino M: The ReWalk powered exoskeleton to restore ambulatory function to individuals with thoracic-level motor-complete spinal cord injury. *Am J Phys Med Rehabil* 2012;00:00Y00.

Objective: The aim of this study was to assess the safety and performance of ReWalk in enabling people with paraplegia due to spinal cord injury to carry out routine ambulatory functions.

Design: This was an open, noncomparative, nonrandomized study of the safety and performance of the ReWalk powered exoskeleton. All 12 subjects have completed the active intervention; three remain in long-term follow-up.

Results: After training, all subjects were able to independently transfer and walk, without human assistance while using the ReWalk, for at least 50 to 100 m continuously, for a period of at least 5 to 10 mins continuously and with velocities ranging from 0.03 to 0.45 m/sec (mean, 0.25 m/sec). Excluding two subjects with considerably reduced walking abilities, average distances and velocities improved significantly. Some subjects reported improvements in pain, bowel and bladder function, and spasticity during the trial. All subjects had strong positive comments regarding the emotional/psychosocial benefits of the use of ReWalk.

Conclusions: ReWalk holds considerable potential as a safe ambulatory powered orthosis for motor-complete thoracic-level spinal cord injury patients. Most subjects achieved a level of walking proficiency close to that needed for limited community ambulation. A high degree of performance variability was observed across individuals. Some of this variability was explained by level of injury, but other factors have not been completely identified. Further development and application of this rehabilitation tool to other diagnoses are expected in the future.

Key Words: Spinal Cord Injury, Paraplegia, Orthotic Devices, Walking, Ambulation, Robotics

The ReWalk Powered Exoskeleton 1

Approximately 200,000 Americans currently live with a disability related to spinal cord injury (SCI).¹ Each year, there are approximately 11,000 new SCIs in this country. Advances in medical care for persons with SCI have resulted in increased rates of survival and longer life span. Many patients with SCI report reduced life satisfaction^{2Y6}; nearly 70% of patients with new SCI have been shown to have reduced quality-of-life, significant levels of pain,⁷ and long-term satisfaction score reduction. Factors associated with reduced life satisfaction among these patients include mobility limitations^{6,8} and changes in occupational roles^{6,8,9} or physical health.^{3,5,6} Some of the major changes contributing to reduced physical health stem from impaired bowel and bladder function, pressure sores, pain, increased muscle tone, obesity, and decreased bone density. Singularly or in combination, these changes can contribute to reduced mobility and social interaction.

Locomotion Options for Individuals with Thoracic-Level Motor-Complete SCI

Manual wheelchair propulsion is the most common mode of locomotion for those with a motor-complete thoracic-level SCI.¹⁰ In wheelchairs, individuals are often able to navigate around a properly modified home and workplace, perform many activities of daily living, and engage in some social and recreational activities. However, because of significant architectural and environmental restrictions, there are limitations associated with wheelchair use.¹¹ Furthermore, the typical seated wheelchair posture does not load the legs in a normal manner¹² and can lead to pressure sores^{13,14} and promote joint contractures. Dependence on arms for propulsion can result in shoulder overuse

syndrome^{15,16} as well as neuropathy. Finally, most wheelchairs do not provide a means for eye-to-eye social interaction with able-bodied adults.

Other modes of locomotion that have been previously used come up far short of even the standard set by the wheelchair. Classic knee-ankle-foot orthoses are heavy, cumbersome, and difficult to don, with most patients discontinuing their use within a short time. Reciprocal gait orthosis, hip guidance, and isocentric reciprocal gait orthoses (IRGO) have brought about improvement over knee-ankle-foot orthoses, allowing increased ambulation for short distances, but long-term adoption is still limited because of high energy consumption from the inefficient gait pattern required and potential for increase upper limb overuse.¹⁷ Other

obstacles to the use of this type of braces include difficulty donning/doffing, obesity, dependence on others for assistance guarding, and mechanical breakdown of the devices. These braces may also interfere with other functional situations, such as moving to/from a couch or car transfers.^{18,19} The American Paraplegia Society has developed clinical practice guidelines for clinicians who work with patients having SCI. The guidelines suggest using a stander at 1 yr after injury for individuals with motor-complete SCI injuries above T1. For injuries between T1 and below, either a stander or orthoses for therapeutic or functional activities are an option.²⁰

Functional electric stimulation systems may enable patients to ambulate for very limited distances.²¹ The first success in functional electric stimulation-assisted gait in paraplegia occurred in the 1970s, with the systems of Kljajic and Bajd in Slovenia²² and Graupe et al. in the United States.²³ Others used implanted electrodes to elicit a gait pattern.²⁴ The technique has many limitations, requiring functioning lower motor neurons for neuromuscular excitability and complete sensory loss to tolerate the significant electrical stimulus needed to achieve an adequate muscular contraction. In addition, the electrical stimulation differs greatly from the physiologic nerve impulse because in functional electric stimulation, all motor units in a muscle group are stimulated simultaneously. This rapidly induces muscle fatigue and results in high energy consumption. To try to overcome these limitations, a hybrid orthosis using swing phase-only functional electric stimulation to assist limb advancement²⁵ and allowing periods of no stimulation in stance phase where electromechanical actuators provide stability is in development.

The walking performances of all the above methods remain quite modest in comparison with normal gait. Speed is very low and energy utilization is high, thus limiting their functional use. The problem of afferent feedback and balance remains unsolved, requiring the use of a walker or crutches for stability, regardless of the system used.

Reduction in upright posture, reduced bone loading relative to functional bipedal walking, diminished practicality, and the potential for a range of unwanted side effects are all factors that have been identified as having a significant negative impact on the overall health of individuals with SCI.^{26,30}

Impact of the Loss of Walking for the SCI Population

Nonrecreational walking accounts for a significant fraction of activity for the average adult. Lower

2 Esquenazi et al. Am. J. Phys. Med. Rehabil. & Vol. 00, No. 00, Month 2012 Copyright

© 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

physical activity levels have been observed after SCI²⁶ resulting from lost motor function, lack of training during acute rehabilitation, decreased access to exercise facilities and adequate adaptive fitness equipment, limited time, and psychologic factors.³¹ The physical deconditioning resulting from largely sedentary lives of individuals with thoracic or higher level SCI is well documented.^{27,29} Lack of standing, ambulation, muscle activity, weight-bearing, and neuroendocrine changes all contribute to rapid and marked alteration in body composition such that muscle mass and bone density decline and adipose tissue increases. Additional complications such as muscle atrophy,³² joint contractures, pressure sores, osteoporosis,³³ increased spasticity, pain, edema, urinary and intestinal stasis, and abnormalities of carbohydrate,³⁴ lipid, and protein metabolism³⁵ may be present.

Exercise has been shown to be an effective contributor to overall health, bone density, a proper level of muscle tone, cardiovascular fitness, regular bowel and bladder function, reduced risk for obesity, heart disease, and reduction in type II diabetes for patients with SCI.^{28,29} Therapeutic exercise for individuals with SCI have several limitations, that is, difficulty in execution,^{36,37} insufficient cardiovascular stimulus, potential for injury, and the need for specialized equipment. Upper body exercises including hand ergometry or weight lifting may increase the likelihood of overuse of already taxed upper limb joints, particularly for thoracic-

level injured individuals who rely on wheelchair propulsion for locomotion and their arms for transfers.^{16,38,39} Having an SCI results in a 40% reduction in work and leisure time compared with preinjury status,⁴⁰ making the adherence to an exercise routine that is separate from, rather than a part of, the activities of daily living more difficult. As an example, in one small study, individuals with SCI who reported relief from pain from a 9-mo exercise program showed a nearly 50% drop in adherence to the program and a resurgence of pain just 3 mos after study conclusion.⁴¹ Functional walking is an excellent means to accomplish exercise without requiring extra time commitments. But this is a difficult option, particularly for those with motor-complete SCI at the thoracic or higher levels.

Description of the ReWalk

The ReWalk is a lower limb powered exoskeleton that allows thoracic or lower level motor-complete individuals with SCI to walk independently (Fig. 1). ReWalk contains independently controlled bilateral hip and knee joint motors, rechargeable batteries, and a computerized control system carried in a backpack. ReWalk users control their walking through subtle trunk motion and changes in center of gravity. A tilt sensor determines the angle of the torso and generates a preset hip and knee displacement

(angle and time) that results in a step. The ankles use a simple double-action orthotic joint with limited motion and spring-assisted dorsiflexion




FIGURE 1 Schematic of the ReWalk I Exoskeleton System. The tilt sensor is used to signal the on-board computer when to take the next step. Body interface is adjustable and uses soft closures. Actuation unit has individual motors for each knee and hip. Wrist unit provides wireless remote access to selected functions in the system. Backpack contains microprocessor, main battery, and backup power supply.

www.ajpmr.com The ReWalk Powered Exoskeleton 3

Copyright © 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

adjustable through screw tension. The device is easily adjustable in height and width and has padded interfaces for calves and thighs and a rigid pelvic frame linking the limbs. Velcro closures with pads, shoes, and a waist belt are used to secure the user in the exoskeleton. Crutches provide standing stability. The subject can interact remotely with the system with a user-operated wrist pad controller that can command sit to stand, stand to sit, and walk activation. The novelty of ReWalk is in the unique manner in which the user is actively involved in controlling walking. The specially designed software algorithm interprets a signal from the torso tilt sensor and generates alternating limb-coordinated motion to produce bipedal walking. At the same time, the system prevents two sequential steps of the same leg. During training, joint angle displacements for the knee and hip can be adjusted using an external computer to optimize the walking characteristics or implement a training mode. A manual mode of operation can

be used to trigger steps bypassing the tilt sensor. The same mode of operation can be used to trigger sit-stand-sit transfers. ReWalk is suitable for adults who have preserved bilateral upper limb function as well as the capacity for assisted standing (such as with a standing frame or braces and crutches). Because the system is battery powered and completely untethered and individuals are fully in control of when they step, ReWalk offers a real option to improve upon the current ambulation standard for individuals with thoracic-level motor-complete SCI. Moreover, because walking in ReWalk is similar to upright bipedal walking, it may offer the potential to overcome some of the physical and psychosocial problems caused by loss of natural walking. In this article, the authors report the findings of their 2-yr clinical trial to assess the safety and efficacy of ReWalk, along with insights into some of its other potential benefits.

METHODS

This was a prospective, single-intervention, open, nonrandomized, partially sponsored by industry trial to evaluate the safety and usability of the ReWalk. The research was approved by the institutional review board at the Albert Einstein Healthcare Network. Informed consent was obtained from all participating subjects before initiating any study-related activity. All subjects had gait training using the ReWalk as the intervention. Subjects were trained for up to 24 sessions of 60- to 90-min duration over approximately 8 wks (target was three times per week). Near the conclusion of the training,

subjects had a performance evaluation visit that consisted of a 6-min walk test, a 10-m walk test, and a gait laboratory evaluation including three-dimensional motion capture and temporospatial data. Dynamic electromyogram was collected during walking from selected lower limb proximal muscles to confirm that no lower limb muscle activity was present. Twelve subjects completed the gait training portion of the study. All subjects who participated in the intervention portion of the trial were followed for approximately 1 yr after active intervention phase to assess any long-term effects; three subjects remain in the follow-up phase.

Subject Selection

Adults with chronic (at least 6 mos postinjury) motor-complete cervical and thoracic (C7-T12) SCI according to American Spinal Injury Association guidelines were targeted for recruitment into the study. Subjects were initially screened by telephone to determine potential eligibility to participate; subsequently, subjects were scheduled to come to the facility to be consented for complete evaluation of their ability to participate. The intake evaluation included dual-energy x-ray absorptiometry, electrocardiography, and leg long bone and lumbar spine x-rays to confirm joint integrity and absence of unhealed fractures or heterotopic ossification that may impede walking. At the time of screening, subjects had to report a history of standing (either with lower limb bracing or a standing frame) on a frequent basis to proceed. Furthermore, a complete neurologic evaluation by a study co-primary investigator (M.S.) was used to confirm injury level; skin integrity; hemodynamic stability; adequate hip, knee, and ankle range of motion; and a spasticity level of 3 or less using the Ashworth scale. The

above information, combined with the absence of osteoporosis on the basis of bone mineral density (BMD) measured from the right-limb femoral neck and the L2 to L4 spine. BMD t scores at the measured sites of greater than -2.5 were required based on the definition of osteoporosis from the World Health Organization.⁴² To make a final determination of subject eligibility to participate, the primary investigator, a board-certified, practicing physiatrist experienced in bio-mechanics of gait and with understanding of the ReWalk technology, reviewed and considered all the available information. Table 1 summarizes inclusion/ exclusion criteria used for this study.

Intervention Phase

After measurement of the subject and fitting with the device, the training program commenced,

4 Esquenazi et al. Am. J. Phys. Med. Rehabil. & Vol. 00, No. 00, Month 2012 Copyright

© 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

TABLE 1 Inclusion and exclusion criteria

Inclusion criteria

Motor-complete cervical (C7Y8) and thoracic

(T1YT12) spinal cord injury/damage according to

American Spinal Injury Association guidelines Male and nonpregnant, nonlactating woman Age 18Y55 yrs

At least 6 mos after injury

Regular use of RGO or KAFOs or able to stand using a standing device (e.g., Easy stand)

Ability to provide informed consent Height of 160Y190 cm

Weight of ≤ 100 kg

Exclusion criteria

History of severe neurologic injuries other than SCI

(MS, CP, ALS, TBI, etc.)

Severe concurrent medical diseases: infections,

circulatory, heart or lung, pressure sores Severe osteoporosis affecting the hip and spine

(t scores ≤ -2)

Severe spasticity (greater than Ashworth 3) or

uncontrolled clonus

Unstable spine or unhealed limb or pelvic fractures Heterotopic ossification that restricts functional

range of motion

Significant contractures (beyond 15 degrees at the

hips or 20 degrees at the knees)
Psychiatric or cognitive comorbidities that may

interfere with the trial

RGO indicates reciprocal gait orthosis; KAFOs, knee- ankle-foot orthoses; SCI, spinal cord injury; MS, multiple sclerosis; CP, cerebral palsy; ALS, amyotrophic lateral sclerosis; TBI, traumatic brain injury.

following a protocol including sit-to-stand and standing activities in the parallel bars, stand-sit transfers, standing balance, and skills related to stepping (Table 2). Training progressed to crutch use placement for balance and advancement. After the acquisition of these underlying skills, the bulk of the training focused on improving and integrating walking performance with step triggering, coordinating step timing and foot clearance, safe and effective stopping, and finally, full self-control using the wrist pad controller. Training sessions were individualized, by the physiotherapist, to facilitate safe and efficient walking. Generally, sessions lasted 75 to 90 mins total, with actual standing and walking training times dependent primarily on equipment adjustment and subject need for rest periods. Heart rate, blood pressure, spasticity, self-reported pain (VAS scale), fatigue (VAS scale), and skin integrity were measured before and after each session by the same member of the research team. Subjects were continually supervised by at least one training staff member during all training. A physical therapist, physical therapist aide, and a biomedical engineer comprised the training staff. Staff received training in the symptomatic detection of autonomic dys-

flexia and orthostatic hypotension and the initial interventions to manage it. Emergency medical staff was available on site to respond at all times during ReWalk use.

Follow-Up

Subjects completed a survey describing their experience at the conclusion of the active intervention phase. This survey included questions about subjects' overall comfort and confidence using ReWalk (Appendix 1) as well as specific inquiries about pain, fatigue, spasticity, intestinal function, and respiratory problems. They were also evaluated at approximately 12 to 15 mos after the intervention to determine whether there were any long-term physical and emotional effects of device use. This consisted of a physical examination, assessment of spasticity and pain, as well as the SF-36v2 Health Survey Questionnaire.

RESULTS

Seventeen subjects were screened, and 13 were considered eligible to participate in the study. Three of those ineligible had severe osteoporosis in the hip and spine and one was too large to comfortably use the ReWalk device. One eligible subject was not able to participate after screening was completed because of study-unrelated personal reasons. Twelve subjects participated in this study. The demographics of these 12 subjects are shown in Table 3. All subjects in the study had thoracic-level motor-complete

TABLE 2 General training schedule Schedule of Gait TrainingVMajor Activities

Session 1: complete necessary measurements Session 1: fitting of device to subject

Session 1: sit to stand between parallel bars (manual

trigger mode)

Sessions 1Y4: putting on device

Sessions 1Y4: sit-to-stand and stand-to-sit within

parallel bars

Sessions 1Y4: stable standing between parallel bars

(without support)

Sessions 1Y4: walking between parallel bars (manual

trigger mode)

Sessions 2Y4: walking between parallel bars (subject

triggered tilt sensor)

Sessions 2Y4: sit-to-stand with help of assistant Sessions 2Y4: basic training in the use of crutches

Sessions 2Y5: walking with crutches (manual trigger

mode)

Sessions 3Y18: walking with crutches using tilt sensor Sessions 3Y18: sit-to-stand and stand-to-sit using

crutches

www.ajpmr.com

The ReWalk Powered Exoskeleton 5

Copyright © 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

TABLE 3 Subject demographics

Age, Weight, Height,

Time Since Injury, yrs

4.8 1.0 2.4 3.3

24.3 13.7 2.1 3.1 19.6 6.0 2.9 6.0

Cause of Injury

Vascular MVA MVA Fall MVA Assault MVA MVA Fall MVA GSW MVA

Injury No. of Level Visits

T8 25 T3 22 T8 26

T11 24 T10 25 T10 24

T4 26

T6 13 T10 23 T12 25

T7 24 T4 26

Subject Code

1

2

3

4

5

6

7

8

9

10

11

12

yrs Sex kg cm

49 Female 68.2 163 45 Male 79.1 183 36 Male 71.4 175 55 Male 81.8 180 46 Male 75.0 185 41
Male 88.6 183 18 Female 56.8 168 36 Male 86.4 185 40 Male 70.5 175 36 Female 56.8 168 20
Male 79.5 188 32 Female 70.5 175

MVA indicates motor vehicle accident; GSW, gunshot wound.

SCI and were nonambulatory without wheelchair before initiating the study. By completion of the trial, all subjects had walked under their own control, without human assistance while using the ReWalk, for at least 50 to 100 m continuously and for a period of at least 5 to 10 mins. Velocities ranged from 0.03 to 0.45 m/sec (mean, 0.25 m/sec) (Table 4). Velocities were considerably slower for the first 2 of the 12 subjects who participated. It is unclear to what extent the authors' own unfamiliarity with the device or not having a fully developed training protocol contributed to the limited outcomes of these first two subjects enrolled in the study. Differences in velocity and distance performance may also be related to level of injury; three groups appear when performance data are stratified by level of injury (Fig. 2).

During the thousands of steps taken through the authors' facility hallways and even outdoors, there was not a single study-related serious event: no falls, no bone fractures, or episodes of autonomic dysreflexia occurred.

Of the 12 subjects, 5 had a definite or possible study-related mild to moderate adverse event. The events consisted of (1) skin abrasions in areas of contact with the device, (2) lightheadedness, and (3) edema of the lower limbs. They were all managed by the appropriate use of (1) foam and padding, (2) caffeinated beverage intake and adjustment of blood pressure medication, (3) elastic stockings and rest, respectively.

There were no detrimental changes in vital signs or complaints of light-headedness with prolonged standing. A summary of changes in spasticity for the bilateral knee and ankle were derived from all sessions for each subject completing the trial (Table 5). The differences were measured from immediately before to immediately after each test session. The scores show

that, in general, spasticity seemed to be favorably impacted by using the ReWalk. Although not formally recorded, total time in standing and walking per session ranged approximately between 20 and 60 mins. Furthermore, subjects 4, 8, and 10 reported in the survey at the end of the training intervention that their overall spasticity improved (noted in Table 5), whereas no subjects indicated that their spasticity had increased.

Pain was measured before and after each test session. Five of the subjects reported a combined 28 times that pain was noticeably reduced; one subject indicated that pain was repeatedly increased immediately after the training intervention. The increased pain complaints were not dissimilar to that produced by standing with a standing frame and, as such, was not considered a ReWalk adverse effect.

TABLE 4 Objective walking performance measures

Subject Code

6-min 6-min Test 10-m Distance, Velocity, Time,

m m/sec min:sec

10-m Test Velocity, m/sec

0.03 0.03 0.25 0.42 0.40 0.45

1

2

3

4

5

6

7

8

9

10 113.5 11 86.7 12 39.6

15.3 0.04 5:10 10.8 0.03 5:25 63.4 0.18 0:40

122.0 0.34 0:24 106.7 0.30 0:25 150.4 0.42 0:22

54.3 0.15 0:55 0.18

aaaa

90.2

^aSubject was unable to participate in this test because of factors unrelated to the study.

0.25 0:39 0.26 0.32 0:31 0.32 0.24 0:36 0.28 0.11 1:09 0.14

6 Esquenazi et al. Am. J. Phys. Med. Rehabil. & Vol. 00, No. 00, Month 2012 Copyright

© 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.




FIGURE 2 Walking performance across injury levels. Figure illustrates the performance on the 6-min and the 10-m walk tests by level of injury and suggests the effect of level of injury on performance (i.e., lower level of injury higher walking velocity and longer distance walked).

In addition to the above data obtained from each test session, self-reported data were collected upon trial termination. The following trends were observed: 3 of 11 subjects reported improved spasticity.

0 of 11 subjects reported that use of ReWalk caused

pain.

1 of 11 subjects reported that use of ReWalk caused

fatigue.

5 of 11 subjects reported improved bowel regulation.

The study therapist and investigators repeatedly and consistently observed all study subjects perspiring sometimes considerably during initial training with the ReWalk. Moreover, cursory review of heart rate data support that elevated heart rates and a modest increase in blood pressure sometimes resulted from use of the ReWalk (Table 5). This improved over time as subjects' fitness and familiarity with using the ReWalk increased.

Finally, a follow-up examination and questionnaire were administered approximately 12 mos after final device use. The general health status, as measured by study clinicians, did not change. Two subjects reported study-unrelated complications in this

phase. One was a victim of stabbing and the other developed a grade 3 pressure sore that required surgery. Evaluation of long-term changes in Modified Ashworth scores, spasm frequency, and interpreted SF-36v2 scores is pending. All subjects had strong positive comments regarding the emotional/ psychosocial benefits of participating in the trial. The authors aim to further evaluate this aspect via the SF-36v2 when long-term follow-up is completed for all study subjects.

Two subjects had implanted medical devices at the time of enrollment: one had a spinal cord stimulation system and another had an intrathecal drug delivery system. Neither patient had device-related complications during the course of the ReWalk intervention. System parameters for both devices were not changed during the course of the training.

DISCUSSION

In a population of thoracic SCI subjects, the ReWalk was found to be safe and effective in allowing independent walking at different levels of ambulatory performance. With careful and regular

TABLE 5 Effect of ReWalk on leg spasticity (Ashworth score), HR, and BP

Subject No. of Times Spasticity No. of Times Spasticity Average Change in BP, Average Change in HR,

Code Improved Got Worse mm Hg beats/min

1 0 0 11/3 4 2 8 18 j3/j2 7 3 40 12 11/2 j5 4 12^a 1 2/j3 7 5 0 0 9/16 19 6 18 33 0/2 6

7 42 8 3^a 9 6 10 1^a 11 0 12 0

0 11/10 26 02/37 1 j3/j2 4 02/04 0 6/4 12 0 6/8 10

^aSubjects for whom self-reported overall spasticity improved. HR indicates heart rate; BP, blood pressure.

www.ajpmr.com

The ReWalk Powered Exoskeleton 7

Copyright © 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

monitoring of walking as described by the manufacturer, use of the ReWalk triggered no serious adverse events and only a relatively small number of minor adverse events. All of these events were straightforward to manage by padding of contact surfaces or physical measures and did not lead to any longer term issues or cause any more than a small delay in trial progress. Subjects did occasionally lose balance and either caught themselves with their crutches or were stabilized by the physical therapist or physical therapist aide. None fell to the ground at any time during walking, standing, or transfers. Subjects were not allowed to use the device unsupervised. Spasticity was not

considered an adverse effect because it is well documented that fatigue may increase this phenomenon and, as such, is a possible expected effect of the training.

With regard to the BMD screening measurements, the community standard for bone densitometry includes the femoral neck and lumbar spine, and these were the locations used in this study. Although it is well recognized that the pattern of bone loss after SCI is different (i.e., the distal femur demonstrating the most profound loss), the authors selected the standard measurement sites because many subjects traveled long distance to participate in the study after obtaining their bone densitometry in their community. The observation that 13 of the 17 screened subjects had nearly normal or only osteopenic levels of BMD at the lumbar spine and femoral neck may not represent the data in the general paraplegic population. The frequent upright positioning required to participate in this study may have attenuated the BMD loss in these subjects. There was no consistent pattern of decreased BMD at either the lumbar spine or the femoral neck in the screened participants. The safety of ReWalk in individuals with osteoporosis has not been formally evaluated at this time, but its design does provide support to the distal femur.

The different levels of performance refers to the range of walking capacity varying from limited household to independent community functional walking. Two of the subjects were particularly limited and unable to sustain 0.1 m/sec walking speed over a 10-m span (Table 5). It is unclear to what extent the authors' own unfamiliarity with the device or lack of a fully developed training protocol contributed to lower performance in these first two study subjects. In addition, several of the subjects were not able to obtain consistent and comfortable use of the remote wrist controller. They were still in control of their walking through self-triggering of steps; however, they were not able to consistently

change between the walking and sitting routines by themselves as they preferred to keep their hands on the crutches. Although many of the subjects were able to perform independently in the hospital setting, they all had limited ability to deal with extraneous factors, such as significant changes in the walking surface, being amidst excessive foot traffic, and the occasional device step misfire (e.g., device not taking a step even when patient commands one).

The authors do not believe that with the prototype device tested and the limited duration of practice allowed in the trial, the subjects would have been able to navigate all terrain situations independently. Additional training seems to be of benefit as improved performance has been noticed in subjects who have continued training beyond the trial. Future implementation of a more formalized user training program based on level of injury may result in higher functional mobility using ReWalk.

Subject performance in the ReWalk was likely influenced by injury level and seemed to fall into broad categories (Fig. 2). Generally, the T2YT4 subjects seemed to perform distinctly differently from those with lower level injuries. The varying levels of performance are reflected in walking speeds that ranged from 0.03 to 0.45 m/sec. Larger studies are needed to determine the validity of this trend. Adequate limb clearance and consistent stepping were the largest obstacles to successful use. Based on this study, having adequate trunk control seems to be a factor for subjects to consistently and re-

peatedly trigger steps to create a fluid gait. This result is not completely unexpected; use of conventional hip-knee-ankle-foot orthoses in patients with paraplegia is typically easier with more caudal levels of thoracic injury.¹⁰ At the completion of this study, training methods and modifications of the device were still in development. The manufacturer has added several changes to the software not available during the trial, including the option of a marching-type gait to facilitate consistent stepping training.

Many subjects were able to walk fairly large distances independently and continuously by the end of the active use phase (9150 m). Although energy consumption was not measured during this safety trial, subjects appeared to work harder than able-bodied individuals walking at the same speeds. Fatigue did not seem to be a limiting factor in performance, nor would it be expected to be the limiting factor in a more natural environment. Subject walking skill level seems to be readily improved with additional training time. The test device required minor but important further development, including, but not limited to, (1) slightly

8 Esquenazi et al. Am. J. Phys. Med. Rehabil. & Vol. 00, No. 00, Month 2012 Copyright

© 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

higher consistency of step trigger performance and (2) improved fail-safe modes if users lose balance and need to reset the device. Other improvements that may be welcomed by the target population include quieter motors and a streamlined appearance. The current commercially available Food and Drug Administration-approved device has implemented many of these suggested improvements.

Both the measured and self-reported data seem to support the premise that consistent ReWalk use has a positive impact on pain, spasticity, and bowel and bladder function. None of the trends could be statistically analyzed because of the pilot nature of this trial. None of the observed benefits were reported to remain once active training stopped. The authors note that medication regimen was documented but not controlled and no changes were made to implanted devices during the trial. The reported improvements cannot be unequivocally attributed to the use of ReWalk; however, it is the authors' expectation that more intense and longer term use of ReWalk may support sustained and statistically significant improvements along these and other health dimensions.

Subjects were positive about the ReWalk experience and there was clear emotional well-being reported from walking. This was particularly evident initially when several subjects were visibly overwhelmed with emotion upon initial walking. The follow-up survey indicated that these effects were not sustained after training discontinuation. However, ReWalk use by an individual in an everyday setting may provide both the emotional and physical health benefits previously described, in addition to improving mobility and the opportunity for increased social interaction. Taken together, these changes may contribute to improved life satisfaction.^{4,43} Neural repair and regeneration approaches, to the degree needed to restore function to patients with complete SCI motor injuries and other neurologic conditions, seem far from being commonplace. The ReWalk and other similar technology may keep individuals with SCI ideally prepared to be candidates for more advanced medical treatments as they are developed by

maintaining physical health, particularly neuromusculoskeletal health below the level of the injury.

The authors continue to explore the potential of ReWalk to restore independent walking function in patients with SCI. Future goals are to expand the use of ReWalk to noncomplete SCI and explore use in other neurologic causes of gait dysfunction such as stroke, brain injury, and cerebral palsy, either as an orthosis or as a gait training intervention.

CONCLUSIONS

The potential for ReWalk to be used as an ambulatory powered orthosis for motor-complete SCI patients is evident from this trial; long-term consistent use is promising as the device is developed for personal home use. Subjects were able to safely participate in training sessions up to three times a week for a period of approximately 2 mos, with no falls or occurrences of autonomic dysreflexia and relatively few minor adverse events related primarily to pressure and irritation from the device or from being upright for prolonged periods of time, for those not used to that. Most subjects were able to achieve a level of walking proficiency that was close to that needed for limited community ambulation in an urban setting. The naturalistic mode of exercise provided by ReWalk may improve some of the health problems reported in this population, thereby possibly lessening the inherent risks associated with their management. Data also suggest that immediate gains may be made in several key areas, including pain management, spasticity, and bowel and bladder function. Furthermore, the use of the ReWalk can offer, in the near future, a greater level of independent upright mobility than currently exists for this population.

ACKNOWLEDGMENTS

Dhan Parker was instrumental in carrying out the subject training; Barbara Hirai and Oliver Woods helped to collect some of the evaluation data in the gait laboratory; and Dr Harry Schwartz contributed to the study design. The authors thank Argo Medica for supplying the ReWalk devices used in this study and their staff for technical support. The authors especially thank the subjects and their families who gave of their time and provided emotional inspiration to complete this study.

REFERENCES

1. CDC. Resources for entertainment education content developers. Available at: <http://www.cdc.gov/health-communication/ToolsTemplates/EntertainmentEd/Tips/SpinalCordInjury.html>. Accessed November 2, 2011
2. Post M, deWitte L, van Asbeck FW, et al: Predictors of health status and life satisfaction in spinal cord injury. *Arch Phys Med Rehabil* 1998;79:395Y401

3. Post M, van Dijk A, van Asbeck F, et al: Life satisfaction of persons with spinal cord injury compared to a population group. *Scand J Rehabil Med* 1998;30:23Y30
4. Putke J, Richards J, Hicken B, et al: Predictors of life satisfaction: A spinal cord injury cohort study. *Arch Phys Med Rehabil* 2002;83:555Y61
5. Schulz R, Decker S: Long-term adjustment to physical disability: The role of social support, perceived

The ReWalk Powered Exoskeleton 9

www.ajpmr.com

Copyright © 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

control, and self-blame. *J Pers Soc Psychol* 1985;48: 1162Y72

6. Fuhrer M, Rintala DH, Hart KA, et al: Relationship of life satisfaction to impairment, disability, and handicap among persons with spinal cord injury living in the community. *Arch Phys Med Rehabil* 1992;73:552Y7
7. Anke AGW, Stenehiem AE, Stanghelle JK: Pain and life quality within 2 years of spinal cord injury. *Paraplegia* 1995;33:555Y9
8. Dijkers M: Correlates of life satisfaction among persons with spinal cord injury. *Arch Phys Med Rehabil* 1999;80:867Y76
9. McColl M, Stirling P, Walker J, et al: Expectations of independence and life satisfaction among ageing spinal cord injured adults. *Disabil Rehabil* 1999;21:231Y40
10. Kirshblum S, Priebe M, Ho C, et al: Spinal cord injury medicine. 3. Rehabilitation phase after acute spinal cord injury. *Arch Phys Med Rehabil* 2007;88:S62Y70
11. Welage N, Liu K: Wheelchair accessibility of public buildings: A review of the literature. *Disabil Rehabil Assist Technol* 2011;6:1Y9
12. Goemaere S, Van Laere M, De Neve P, et al: Bone mineral status in paraplegic patients who do or do not perform standing. *Osteoporos Int* 1994;4:138Y43
13. Wall J: Preventing pressure sores among wheelchair users. *Prof Nurse* 2000;15:321Y4

14. Regan M, Teasell R, Wolfe D, et al: A systematic review of therapeutic interventions for pressure ulcers after spinal cord injury. *Arch Phys Med Rehabil* 2009;90:213Y31
15. Groah S, Lanig I: Neuromusculoskeletal syndromes in wheelchair athletes. *Semin Neurol* 2000;20:201Y8
16. Alm M, Saraste H, Norrbrink C: Shoulder pain in persons with thoracic spinal cord injury: Prevalence and characteristics. *J Rehabil Med* 2008;40:277Y83
17. Chafetz R, Johnston T, Calhoun C: Lower limb orthoses for persons with spinal cord injury, in Hsu J, Michael J, Fisk J (eds): *AAOS Atlas of Orthoses and Assistive Devices*, chapter 27. Philadelphia, PA, Mosby, Elsevier Inc, 2008;357Y72
18. Kaplan L, Grynbaum B, Rusk H, et al: A reappraisal of braces and other mechanical aids in patients with spinal cord dysfunction: Results of a follow-up study. *Arch Phys Med Rehabil* 1966;47:393Y405
19. Franceschini M, Baratta S, Zampolini M, et al: Reciprocating gait orthoses: A multicenter study of their use by spinal cord injured patients. *Arch Phys Med Rehabil* 1997;78:582Y6
20. Consortium for Spinal Cord Medicine: Outcomes following traumatic spinal cord injury: Clinical practice guidelines for health-care professionals. *J Spinal Cord Med* 2000;23:289Y316
21. Stein R, Chong S, James K, et al: Electrical stimulation for therapy and mobility after spinal cord injury. *Prog Brain Res* 2002;137:27Y34
22. Kljajic M, Bajd T, Stanic U: Quantitative gait evaluation of hemiplegic patients using electrical stimulation orthoses. *IEEE Trans Biomed Eng* 1975;22:438Y41
23. Graupe D, Kohn K, Basseas S: Above- and below-lesion EMG pattern mapping for controlling electrical stimulation of paraplegics to facilitate unbraced walker-assisted walking. *J Biomed Eng* 1988;10:305Y11
24. Thoma H, Frey M, Gruber H, et al: First implantation of a 16-channel electric stimulation device in the human body. *Trans Am Soc Artif Intern Organs* 1983;29:301Y6
25. Audu M, To C, Kobetic R, et al: Gait evaluation of a novel hip constraint orthosis with implication for walking in paraplegia. *IEEE Trans Neural Syst Rehabil Eng* 2010;18:610Y8
26. Buchholz A, McGillivray C, Pencharz P: Physical activity levels are low in free-living adults with chronic paraplegia. *Obes Res* 2003;11:563Y70

27. Brenes G, Dearwater S, Shapera R, et al: High density lipoprotein cholesterol concentrations in physically active and sedentary spinal cord injured patients. Arch Phys Med Rehabil 1986;67:445Y50
28. Martin-Ginis KA, Latimer AE, McKechnie K, et al: Using exercise to enhance subjective well-being among people with spinal cord injury: The mediating influences of stress and pain. Rehabil Psychol 2003;48:157Y64
29. Jacobs P, Nash M: Exercise recommendations for individuals with spinal cord injury. Sports Med 2004;34: 727Y51
30. Stevens S, Caputo J, Fuller D, et al: Physical activity and quality of life in adults with spinal cord injury. J Spinal Cord Med 2008;31:373Y8
31. Levins S, Redenbach D, Dyck I: Individual and societal influences on participation in physical activity following spinal cord injury: A qualitative study. Phys Ther 2004;84:496Y510
32. Castro M, Apple DJ, Staron R, et al: Influence of complete spinal cord injury on skeletal muscle within 6 months of injury. J Appl Physiol 1999;86:350Y8
33. Garland D, Stewart C, Adkins R, et al: Osteoporosis after spinal cord injury. J Orthop Res 1992;10:371Y8
34. Elder C, Apple D, Bickel C, et al: Intramuscular fat and glucose tolerance after spinal cord injury: A cross-sectional study. Spinal Cord 2004;42:711Y6
35. Dallmeijer A, van der Woude L, van Kamp G, et al: Changes in lipid, lipoprotein and apolipoprotein profiles in persons with spinal cord injuries during the first 2 years post-injury. Spinal Cord 1999;37:96Y102
36. Thrasher A, Graham G, Popovic M: Reducing muscle fatigue due to functional electrical stimulation using random modulation of stimulation parameters. Artif Organ 2005;29:453Y8
37. Estigoni E, Fornusek C, Smith RM, et al: Evoked EMG and muscle fatigue during isokinetic FES cycling in individuals with SCI. Neuromodulation 2011;14:349Y55
38. Sie I, Waters R, Adkins R, et al: Upper extremity pain in the postrehabilitation spinal cord injured patient. Arch Phys Med Rehabil 1992;73:44Y8
39. Jain N, Higgins L, Katz J, et al: Association of shoulder pain with the use of mobility devices in persons with chronic spinal cord injury. PM&R 2010;2:896Y900

10 Esquenazi et al. Am. J. Phys. Med. Rehabil. & Vol. 00, No. 00, Month 2012

Copyright © 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

40. Schoenherr MC, Groothoff JW, Mulder GA, et al: Participation and satisfaction after spinal cord injury: Results of a vocational and leisure outcome study. Spinal Cord 2005;43:241Y8
41. Ditor DS, Latimer AE, Ginis KAM, et al: Maintenance of exercise participation in individuals with spinal cord injury: Effects on quality of life, stress and pain. Spinal Cord 2003;41:446Y50

42. Kanis J, Melton LI, Christiansen C, et al: The diagnosis of osteoporosis. J Bone Miner Res 1994;9:1137Y41

43. Wood-Dauphine´e S, Exner G, Bostanci B, et al; SCI Consensus Group: Quality of life in patients with

spinal cord injuryVBasic issues, assessment, and- recommendations. Restor Neurol Neurosci 2002;20: 135Y49

APPENDIX 1

Participant Satisfaction Questionnaire

The 10 questions below were asked of each subject upon the completion of the active participation phase of the research project. Subjects were given the following five answer choices to each question:

-
1. The ReWalk training process is easy.
 2. Fitting and adjusting ReWalk is easy.
 3. I felt comfortable using ReWalk.
 4. Using ReWalk did not cause me any pain.
 5. I didn't feel especially fatigued when using ReWalk.
 6. When I finished training, I felt comfortable using ReWalk.
 7. Use of the device has lessened the spasticity in my lower limbs.
 8. I didn't experience any respiratory problems when using the device.
 9. During the period when I was using the device I felt an improvement

in my intestinal activity.

10. At the end of the training period I felt confident using the device.

Strongly Disagree (1)

Disagree (2)

Neither Agree or Disagree (3)

Agree (4)

Strongly Agree (5)

www.ajpmr.com

The ReWalk Powered Exoskeleton 11

Copyright © 2012 Lippincott Williams & Wilkins. Unauthorized reproduction of this article is prohibited.

Amaç:

Bu çalışmanın amacı, omurilik yaralanmasına bağlı paraplejisi olan kişilerin rutin ayakta fonksiyonel işlevleri yerine getirmelerinde ReWalk'un emniyet ve performansını değerlendirmektir.

Tasarı:

Bu, ReWalk ile çalışan dış iskeletin emniyetini ve performansını açıkça, karşılaştırılmayan, randomize olmayan bir çalışmaydı. 12 denegin hepsi aktif müdahaleyi tamamladı; Üçü uzun vadeli izlemde kaldı.

Sonuçlar:

Eğitim sonrasında, tüm denekler, ReWalk'ı kullanırken insan yardımına gerek kalmaksızın, sürekli olarak en az 50 ila 100 m boyunca, en az 5 ila 10 dakika sürekli ve hızları 0.03 ila 0,45 m / sn (ortalama, 0,25 m / sn). Oldukça düşük yürüme özelliklerine sahip iki denek hariç tutulduğunda, ortalama mesafeler ve hızlar önemli ölçüde iyileşmiştir. Bazı kişiler deneme sırasında ağrı, bağırsak ve mesane fonksiyonlarında düzelme ve spastisiteyi bildirdiler. Tüm denekler, ReWalk kullanımının duygusal / psikososyal faydaları konusunda güçlü pozitif yorumlarda bulundu.

Sonuçlar:

ReWalk, motor tam torasik düzeydeki omurilik hasarı hastaları için güvenli bir ambulator ortez olarak önemli bir potansiyele sahiptir. Çoğu konu, sınırlı topluluk ambulasyonu için gereken seviyeye yakın düzeyde bir düzeyde yürüme yeterliliği elde etti. Bireyler arasında yüksek düzeyde bir performans değişkenliği gözlemlendi. Bu değişkenliği bazıları yaralanma seviyesine göre açıklanmıştır, ancak diğer faktörler tamamen tanımlanmamıştır. Gelecekte bu rehabilitasyon aracının diğer teşhislere daha da geliştirilmesi ve uygulanması beklenmektedir.